

# Using Evzio<sup>®</sup> auto-injector

- Use EVZIO<sup>®</sup> exactly as prescribed by your healthcare provider and as described in the *Patient Instructions for Use*.
- Each EVZIO<sup>®</sup> auto-injector contains only one dose of medicine.
- As part of your opioid emergency plan, practice using the Trainer for EVZIO<sup>®</sup>.
- Make sure your family and caregiver also practice using the Trainer for EVZIO<sup>®</sup>.

FIND OUT MORE AT

→ [BEHAVIORALHEALTH.ND.GOV/OPIOIDS](https://behavioralhealth.nd.gov/opioids)

NORTH  
**Dakota**  
Be Legendary.™

Behavioral Health  
HUMAN SERVICES

**OPIOIDS**

FILL with CARE

# Administration Steps

---

- 1 Pull EVZIO® from the outer case.
- 2 Pull off the red safety guard.
- 3 Place the black end against the middle of the patient's outer thigh, through clothing (pants, jeans, etc.) if necessary, then press firmly and hold in place for 5 seconds.
- 4 After using EVZIO®, get emergency medical help immediately.

FOR MORE INFORMATION, GO TO  
[EVZIO.COM/PATIENT/ABOUT-EVZIO](http://EVZIO.COM/PATIENT/ABOUT-EVZIO)